


# PEDIATRIC ANAL DILATORS - SINGLE PATIENT USE



## SPECIFIC

- **Smooth finish** for patient comfort.
- **Heat stable** to 270° F, may be boiled or cold sterilized.
- **Size imprinted** on each dilator.
- **Finger grip** for more secure handling.
- **Wide size range**, from 3mm through 19mm in 1mm increments.
- **Available individually** or in sets of eight in any size(s) of your choice.
- **Re-closeable packaging** for convenient storage and patient distribution.
- Not made with natural rubber latex.
- **Custom sizes** may be available upon request.

