


# PEDIATRIC ANAL DILATORS - SINGLE PATIENT USE



## SPECIFIC DESIGN POINTS:

- **Smooth finish** for patient comfort.
- **Heat stable** to 270° F, may be boiled or cold sterilized.
- **Size imprinted** on each dilator.
- **Finger grip** for more secure handling.
- **Wide size** range, from 3mm through 19mm in 1mm increments.
- **Available individually** or in sets of eight in any size(s) of your choice.
- **Re-closeable packaging** for convenient storage and patient distribution.
- **Latex free.**
- **Custom sizes** may be available upon request.


# PEDIATRIC ANAL DILATORS

## ORDERING INFORMATION

## SUPPLIED NON-STERILE

CATALOG #	DESCRIPTION
-----------	-------------

PAD03	Pediatric Anal Dilator, 3mm
PAD04	Pediatric Anal Dilator, 4mm
PAD05	Pediatric Anal Dilator, 5mm
PAD06	Pediatric Anal Dilator, 6mm
PAD07	Pediatric Anal Dilator, 7mm
PAD08	Pediatric Anal Dilator, 8mm
PAD09	Pediatric Anal Dilator, 9mm
PAD10	Pediatric Anal Dilator, 10mm
PAD11	Pediatric Anal Dilator, 11mm
PAD12	Pediatric Anal Dilator, 12mm
PAD13	Pediatric Anal Dilator, 13mm
PAD14	Pediatric Anal Dilator, 14mm
PAD15	Pediatric Anal Dilator, 15mm
PAD15L	Pediatric Anal Dilator, 15mm Long
PAD16	Pediatric Anal Dilator, 16mm
PAD16L	Pediatric Anal Dilator, 16mm Long
PAD17	Pediatric Anal Dilator, 17mm
PAD17L	Pediatric Anal Dilator, 17mm Long
PAD18	Pediatric Anal Dilator, 18mm
PAD18L	Pediatric Anal Dilator, 18mm Long
PAD19	Pediatric Anal Dilator, 19mm
PAD19L	Pediatric Anal Dilator, 19mm Long
PADS1	**Pediatric Anal Dilator Set, 7mm-14mm



Order for Dr.: \_\_\_\_\_

Patient: \_\_\_\_\_

Required By: \_\_\_\_\_

**Catalog #:**

**Qty:**

_____	_____
_____	_____
_____	_____
_____	_____
_____	_____
_____	_____
_____	_____
_____	_____
_____	_____
_____	_____

\*\*Set may be custom made; select eight sizes of your choice.

\$9.00/each when purchased in quantities of eight or more (\$72.00/set).  
 \$10.50/each when purchased in quantities of fewer than eight.

**Limited Warranty, Limitation of Liability and Disclaimer of Other Warranties:**

SSP warrants that reasonable care was used in the manufacture and production of this product. Because SSP has no control over the conditions of use, patient selection, surgical procedure, post-surgical stresses, or handling of the device after it leaves our possession, SSP does not warrant either a good effect or against an ill effect following

its use. SSP shall not be responsible for any incidental or consequential loss, damage or expenses directly or indirectly arising from the use of this product. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to any implied warranties of merchantability or fitness for use.



**Toll Free 888.878.0811**

406.961.0102 · Fax: 406.961.0103 · 1131 North U.S. Hwy 93 · Victor, MT 59875  
 email: info@ssp-inc.com